

Brindlewythe Biosciences Ltd

OB11 antagonist backup set: five-CYP panel

Final study report | Revision B

Brindlewythe Biosciences Ltd

Study BW19_0436

Accession B-ACC-2019-07416

Client: Orvessant Biosciences Pty Ltd

Date: 29 Nov 2019

Author: Dr B. Rowcroft, Screening Operations Director

1 Summary and plate acceptance

40 compounds were tested at 10 μ M against CYP1A2, CYP2C9, CYP2C19, CYP2D6 and CYP3A4 using human liver microsomes and probe substrates.

Reference inhibitor responses supported plate acceptance for all probe systems. Measured control IC₅₀ values are shown with historical comparators. Historical values are reference points, not a requirement for numerical identity.

Table 1. Study benchmark

Metric	Value
Project	BW19_0436
Client	Orvessant Biosciences Pty Ltd
Compounds tested	40
Series / chemotype	Orvessant piperidyl benzamide series 11 (Piperidyl benzamide)
Assays	CYP inhibition panel (1A2/2C9/2C19/2D6/3A4) at 10 μ M
Turnaround (calendar days)	17
Price	GBP 11,500
CYP1A2 >50% inhibition @10 μ M	7
CYP2C9 >50% inhibition @10 μ M	10
CYP2C19 >50% inhibition @10 μ M	9
CYP2D6 >50% inhibition @10 μ M	13
CYP3A4 >50% inhibition @10 μ M	9
Compounds with any CYP >50%	17
Ketoconazole (3A4) IC ₅₀ (μ M)	0.048

Control acceptance precedes interpretation. The benchmark summarises the study; the complete compound tables below are the basis of selection. Screening decisions are recorded separately from the endpoints.

2 Principal results

Table 2. Complete screening set: percent inhibition at 10 μ M

Compound	CYP1A2 (%)	CYP2C9 (%)	CYP2C19 (%)	CYP2D6 (%)	CYP3A4 (%)
OB11-08900	2	22	-4	81	4
OB11-08904	4	16	7	-4	8
OB11-08906	23	80	20	22	56
OB11-08912	15	3	-5	28	2
OB11-08916	-11	10	25	59	42
OB11-08918	70	100	100	3	84
OB11-08924	7	16	4	11	20
OB11-08928	1	1	31	15	0
OB11-08930	15	4	8	2	9
OB11-08934	-2	11	20	-2	4
OB11-08936	31	30	74	68	65
OB11-08940	16	17	22	70	7
OB11-08946	31	62	30	77	35
OB11-08954	8	37	44	13	3
OB11-08960	100	90	100	98	93
OB11-08978	9	3	22	92	2
OB11-08982	6	3	-1	2	5
OB11-08984	9	6	1	3	8
OB11-08990	72	95	100	92	88
OB11-08994	7	7	-2	15	0
OB11-09000	18	5	2	14	1
OB11-09006	-0	8	3	50	28
OB11-09012	4	17	32	-3	14
OB11-09014	1	10	36	-3	3
OB11-09020	56	86	32	63	49
OB11-09024	7	10	25	-3	14
OB11-09026	20	6	41	4	3
OB11-09032	96	98	100	40	98
OB11-09036	12	15	50	-0	2
OB11-09038	16	32	48	20	19
OB11-09042	0	1	4	9	-1
OB11-09048	3	-3	36	14	4
OB11-09060	9	18	23	96	35
OB11-09068	13	30	48	84	-9
OB11-09072	27	92	92	1	85
OB11-09074	13	4	4	4	2
OB11-09078	-4	-5	2	1	1
OB11-09080	93	98	93	56	91
OB11-09086	97	92	78	79	92
OB11-09090	-4	10	12	24	24

Single concentration 10 μ M, duplicate incubations, mean reported; negative values = no inhibition

Controls (IC₅₀, μ M): furafylline (CYP1A2) 1.8; sulfaphenazole (CYP2C9) 0.35; ticlopidine (CYP2C19) 0.88; quinidine (CYP2D6) 0.050; ketoconazole (CYP3A4) 0.048

17 compounds exceeded 50% inhibition for at least one isoform. CYP2D6 had the largest number of flags. Isoform counts overlap and must not be added to obtain a compound count.

Compound identifiers follow the sample cards. Units are part of the headings. Compare like endpoints under the stated conditions. A position in a table is not an estimate of clinical performance.

3 Study remit and chronology

Orvessant Biosciences Pty Ltd commissioned 5-isoform CYP panel to support a chemistry screen gate for this screening set. The remit was not a complete disposition or safety assessment.

Study work began on 12 Nov 2019. The interim checkpoint was 20 Nov 2019. This final revision closes the study on 29 Nov 2019. Dr B. Rowcroft is responsible for scientific release; Ines Alvarado checked plate acceptance and table transfer.

The completed study covers 40 compounds under accession B-ACC-2019-07416. The accession links submitted material, sample cards and run sheets. No later screening set is included.

4 Test-lot handling

The sample-card identifier links each test lot across all probe systems. A shared identifier permits comparison of the panel pattern without implying that different probe reactions are interchangeable. Handling observations remain with the run sheet.

Sample-card identifiers carried through preparation, assay allocation and result collation. Working solutions were prepared for the scheduled assay with matched vehicle. Results refer to the submitted test lot, rather than a claim that a later resynthesis will behave identically.

Preparation details, positions and handling observations remain in the run sheets. Qualifications stay with the endpoint; they are not removed to simplify a ranking.

5 Assay procedure

Human liver microsomes were incubated with an isoform-specific probe substrate and test compound at 10 μ M. Vehicle wells supplied the uninhibited response. Reference inhibitors were run in concentration response for the corresponding probe systems. Vehicle content was matched across test and control wells.

Reactions were initiated with the cofactor system and terminated for measurement of probe product. Incubation conditions were held constant within each isoform assay. Percent inhibition was calculated relative to the matched vehicle response after background correction.

The panel used the same nominal test concentration throughout. It provides a screening interaction map, rather than a test-compound IC50 determination. No preincubation comparison designed to assess time-dependent inhibition was included.

6 Control assessment and release

Reference inhibitor responses supported plate acceptance for all probe systems. Measured control IC₅₀ values are shown with historical comparators. Historical values are reference points, not a requirement for numerical identity.

Table 3. Reference inhibitor responses

Isoform	Reference IC ₅₀ (μM)	Historical comparator (μM)
CYP1A2	1.8	1.8
CYP2C9	0.35	0.35
CYP2C19	0.88	1.2
CYP2D6	0.050	0.060
CYP3A4	0.048	0.030

Table 4. Screen flags by isoform

Isoform	Compounds above 50% at 10 μM
CYP1A2	7
CYP2C9	10
CYP2C19	9
CYP2D6	13
CYP3A4	9

Control potency is expressed in μM. Compound results are percent inhibition at a stated concentration. These endpoints are not interchangeable; a reference concentration response does not supply an IC₅₀ for a test compound.

Ines Alvarado checked that control assessment preceded the screening interpretation. Dr B. Rowcroft authorised release. Acceptance supports use within the stated assay; it does not eliminate compound-specific questions that need follow-up.

7 Data reduction

Percent inhibition is calculated from probe product relative to vehicle. Flags use strictly above 50% at 10 μ M. Negative reported inhibition is retained as variation around the vehicle response; it is not evidence of enzyme activation. A result below the threshold does not establish absence of inhibition at another concentration.

Reported precision and units are retained. Rounding must not manufacture distinctions between near neighbours. The screening decision should use complete rows with qualifications visible.

Summary endpoints and individual observations have different purposes. The benchmark does not replace the detailed table. No additional fitted values or composite scores have been inferred for presentation.

8 Screening-set interpretation

The ether-linked backup set offered reduced CYP liability in the client's chemistry assessment. This supported exposure comparison, but did not establish an exposure margin or an inhibition mechanism.

A broad interaction pattern and an isolated isoform response raise different follow-up questions. The number of flags is a screening aid, not a quantitative ranking of clinical interaction risk. Read the full compound row before assigning a follow-up.

A chemistry association can identify the next experiment without establishing the molecular cause of the assay response. This report makes no target-mechanism claim from screening data alone.

9 Screen gate

Orvessant commissioned an exposure comparison after the ether-linked backup set reduced the CYP liability enough for its screen gate.

Carry the interaction map into the exposure comparison. Preserve sample-card identity for the selected backups and assess exposure alongside the relevant isoform responses. No exposure result is incorporated into this report.

The screen gate records the client's use of the data. A proposed follow-up is not a completed result. No endpoint from that later work is incorporated into this signed report.

10 Study limitations

A single concentration cannot distinguish steep from shallow concentration responses. Nominal incubation concentration is not a measured free concentration. Binding, solubility and probe substrate choice may affect the apparent response.

The panel does not distinguish reversible, irreversible or time-dependent inhibition. It provides neither induction measurements nor metabolite identification. Clinical interaction margins require additional evidence.

These endpoints apply to the submitted screening set and test lots. No extrapolation is made to an untested analogue, formulation or biological system. Later measurements should retain identifiers and conditions for comparison.

11 Conclusions

5-isoform CYP panel is complete for BW19_0436. Control assessment supports release of the detailed tables for the specified client decision.

Orvessant commissioned an exposure comparison after the ether-linked backup set reduced the CYP liability enough for its screen gate.

The results support the stated next step. They do not alone resolve exposure, efficacy or nomination suitability. Any later integrated assessment must identify its additional supporting studies.

12 Record index and authorisation

The retained study record comprises sample cards linked to B-ACC-2019-07416, preparation and assay run sheets, control assessment, result tables and this final report. The interim checkpoint forms part of the chronology and does not replace the completed dataset.

Revision B is the final report issued to Orvessant Biosciences Pty Ltd on 29 Nov 2019. Read principal results with the methods and limitations. Preserve study identity, units and qualifications when using tables for the screen gate.

Scientific release: Dr B. Rowcroft, Screening Operations Director. Plate acceptance and report-table transfer: Ines Alvarado, Assay Quality Coordinator. Signatures apply to this completed study and its bounded interpretation.

Name	Role	Signature	Date
Dr B. Rowcroft	Screening Operations Director — scientific release	<u>D Rowcroft</u>	29 Nov 2019
Ines Alvarado	Assay Quality Coordinator — plate acceptance and table transfer	<u>I Alvarado</u>	29 Nov 2019